

May 9, 2019

To Whom It May Concern,

Penicillin-Streptomycin solutions are typically a colorless solution. However, this material can range in color from colorless to pink to yellow. This is a normal characteristic and this will not affect the potency or application of the product. The coloring is a carry over from the manufacturing process of Streptomycin. The genus that this is isolated from (Actinomycetes Streptomyces griseus) is responsible for a wide variety of pigments.

Please contact us if you have any questions about the information provided.

Sincerely,



Ling Gao, Ph.D.  
Technical Application Scientist III  
1-800-955-6288  
cellculturesupport@thermofisher.com